Wasser & Reiniger

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Water & Purification

**APPLICATION IN DETAIL**

## Sterile Tankbelüftung

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## Sterile tank venting

### Sterile and water-resistant for safe pressure compensation


*In the fields of pharmaceuticals and food as well as cosmetics, sterile venting and bleeding of the storage tanks is absolutely vital. As tank venting filters, highly porous hydrophobic PTFE membrane filter cartridges guarantee safe pressure compensation during the filling and draining processes, which protects the products in the tank from contaminants and bacteria.*

*Concentrated fruit juices, liquid chocolate, sauces or intermediate products for manufacturing toothpaste must all be stored temporarily in sealed tanks for further processing. In addition to plastics, different thicknesses of stainless steel are used for the walls of tanks. These parameters in addition to the height and diameter of the tank determine the low-pressure resistance of a storage tank to elastic or plastic denting as well as its high-pressure resistance. A ventilation connection is attached to the tank for pressure compensation. In sensitive areas like food, drinks, pharmaceuticals, cosmetics or ultra-pure water, venting and bleeding must be carried out on a sterile basis such that no contaminants or bacteria can enter the storage tank. Sterile tank venting filters are used for this.*

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*Fig. 1: The job of the sterile tank venting system is to allow air to flow in and out without the pressure resistance of the filter creating under- or over-pressure in the tank and without an undesirable exchange of contaminants, germs and bacteria taking place.*

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### Operating principle of the filters

*When filling a storage tank, the supply pump pumps a certain volume of liquid per unit of time into the tank. If there was no venting at the highest point in the tank, the gas or gas/air mixture that it contains would be compressed and, in a worst case scenario, the pressure in the tank could increase leading to irreversible deformation of the tank or even to its bursting. This pressure increase can be explained scientifically using the equation of state of a hypothetical ideal gas ( $p \cdot V = m \cdot R \cdot T$ ) that is also referred to as the general ideal gas law of thermodynamics. In many processes and applications, air is undesirable in storage tanks due to the oxidation reactions of the oxygen that it contains. The liquid is superimposed with a protective gas. Here, the described processes apply in a similar way. To simplify things, only air is taken into consideration in the process.*

*When the product is pumped out of the storage tank, a volume of liquid per unit of time is removed from the tank. According to the general ideal gas law, if no air flowed back into the tank, a vacuum would be created in the tank. In a worst case situation, this would lead to the tank contracting and being irreversibly deformed. This means that the job of the sterile tank venting filter is to allow air to flow in and out without the pressure resistance of the filter creating under- or over-pressure in the tank and without an undesirable exchange of contaminants, germs and bacteria taking place.*

*Since physical reasons dictate that the air pocket in a tank always forms above the liquid, the sterile tank venting facility is always located at the top on the storage tank or at least above the maximum possible fill level. The sterile tank venting facility has connections for the nozzle on the storage tank and a support for the filter element. A protective hood (cartridge filter housing) or a dome protect the filter element from contamination from above. (Fig. 2: On next page)*

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## The problem of condensation water

*Often, the liquid to be stored comes out of the process cooled and is filled directly into the storage tank. In the summer, the tank and the ambient air are hot. The difference in temperature results in condensation on the outside of the tank wall. Water droplets form on the surface of the tank and on the outside of the sterile tank venting filter.*

*On a physical basis, the same thing happens if warm or hot liquid is pumped into the storage tank. In this case, the temperature difference between the liquid and the ambient air and the storage tank leads to condensation forming inside the tank and on the inside of the sterile tank venting filter. In both cases, the condensation liquid can have a seriously negative effect on the sterile tank venting filter.*

*If, in addition, the condensate penetrates the pores of the filter element, it is no longer possible for air to flow through the filter, since its pores are blocked by the water that has penetrated. This means that the function of the sterile tank venting filter would be ineffective. It would no longer be possible for pressure compensation to take place, which could have possibly fatal consequences.*

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*Fig. 2:* Screwing the cover off the cartridge filter housing.

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
*Putting the cover back on and screwing it tight.*

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## Hydrophobic membrane materials

To ensure that condensates cannot negatively affect through flowing of the sterile tank venting filter, the venting filter must be fitted with hydrophobic (water-repellent) filter materials.

Membranes made of polytetrafluoroethylene (PTFE) are particularly suitable for this. PTFE is an extremely hydrophobic material. Due to their water-repellent but vapour-permeable material properties, yarns made of PTFE fibres are widely used as fabrics for outdoor and ski jackets as well as for tents. This material property is exploited with sterile tank venting filters. Condensate droplets stay on the surface of the membranes and do not enter the pores of the PTFE membranes. Rather, they drip off and return to the tank or are drained via an external condensate drain on the storage tank.

This means that a dripping sterile tank venting facility is a good sign! (Fig. 4 Illustration of Sterile tank venting)

In conjunction with WFPTFE hydrophobic membrane filter elements, Wolftechnik WTBE housings made of polypropylene (PP) or stainless steel are (Fig. 3) used for dust-free and sterile venting and bleeding of storage tanks. The two-piece design consists of a screw-off hood and the support with a tank nozzle connection as a DN40 GF screw fitting. The housings are available in three sizes to hold a WFPTFE membrane filter element with adapter configuration Code 0 in lengths of 5", 10" or 20", which means that they are suitable for throughputs of 5-50 Nm<sup>3</sup>/h at differential pressure of 10-25 mbar.

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**Fig. 3: Hydrophobic membrane filter elements**

PTFE membranes with a filter fineness of  $\leq 0.2\mu\text{m}$  safely hold back bacteria and germs. Due to the hydrophobic properties of the highly porous membranes, the tank venting filters guarantee safe pressure compensation.

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*PTFE membranes with a filter fineness of  $\leq 0.2\mu\text{m}$  safely hold back bacteria and germs. WFPTFE membrane filter elements with hydrophobic PTFE membranes can be supplied with a filter fineness of  $0.05\mu\text{m}$  to  $1.0\mu\text{m}$ . The membrane material has a highly porous structure and homogeneous pore distribution. All of the components of the support structure, adapter and end caps are made of polypropylene and are joined thermally to the PTFE membranes in a clean room without using additives or surface-active substances. The filter surface is  $0.9\text{ m}^2$  per  $10''$  element which is well above the standard of similar products. This ensures a long service life (long life cycle, high contaminant absorption capacity) and low differential pressures even at high throughput capacities.*

### Ratings of filter housings and filter cartridges

*The rating of sterile tank venting filters based on PTFE membrane filter elements with a filter fineness of  $0.2\mu\text{m}$  and market standard gas throughput capacities complies with the maximum permissible differential pressure for the storage tank (pay attention to the manufacturer's instructions!) In general, this is 10mbar or 25mbar. Relative to approximately  $1\text{ m}^2$  of PTFE membrane filter material of  $0.2\mu\text{m}$ , a gas throughflow of approximately  $10\text{ Nm}^3/\text{h}$  at 10mbar or  $25\text{ Nm}^3/\text{h}$  at 25 mbar is to be expected (see illustration below)*

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